## Tutoriel Phyphox et Arduino

Frédéric Bouquet  
Gautier Creutzer

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5

```
#include <phyphoxBle.h>
#include <Arduino_LSM9DS1.h>
#include <Arduino_LPS22HB.h>
#include <Arduino_HTS221.h>
#include <Arduino_APDS9960.h>

char board_name[] = "Phyphox1"; // to change the name displayed by the board using BLE: nano_arduino_33_ble
float choice = 0.0;
float accx, accy, accz, acc;
--
```

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#### Arduino Nano 33 BLE Sense

AAc	Absolute Acceleration Plots absolute acceleration measure...
AGy	Absolute Angular Velocity... Plots absolute angular velocity mea...
Acc	Acceleration Plots acceleration measured by the...
Gyr	Angular Velocity (Gyroscop... Plots angular velocity measured by t...
Ill	Illuminance Plots illuminance measured by the A...
Mag	Magnetic Field Plots magnetic field measured by th...
Pre	Pressure Plots pressure measured by the Ard...
Tem	Temperature and Humidity Plots temperature and humidity mea...

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Acc	Acceleration Plots acceleration measured by the Arduin...
Pre	Analog Input Plots A0 analog input.
Gyr	Angular Velocity (Gyroscope) Plots angular velocity measured by the Ard...
Ill	Illuminance Plots illuminance measured by the Arduino...
Mag	Magnetic Field Plots magnetic field measured by the Ardu...
Pre	Pressure Plots pressure measured by the Arduino c...
Tem	Temperature and Humidity Plots temperature and humidity measured...

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