

# Phyphox and Arduino Tutorial

Frédéric Bouquet  
Gautier Creutzer

## Physics Reimagined

Laboratory of Solid State Physics  
Paris-Saclay University

October 2020

## Download the Arduino IDE



### ARDUINO 1.8.13

The open-source Arduino Software (IDE) makes it easy to write code and upload it to the board. It runs on Windows, Mac OS X, and Linux. The environment is written in Java and based on Processing and other open-source software.  
This software can be used with any Arduino board. Refer to the [Getting Started](#) page for Installation instructions.

**Windows** Installer, for Windows 7 and up  
**Windows** ZIP file for non admin install

**Windows app** Requires Win 8.1 or 10  
[Get](#)

**Mac OS X** 10.10 or newer

**Linux** 32 bits  
**Linux** 64 bits  
**Linux** ARM 32 bits  
**Linux** ARM 64 bits


[Release Notes](#)  
[Source Code](#)  
[Checksums \(sha512\)](#)

Download and install the Arduino IDE on [their website](#).

Download our program on [our site](#).


Connect the board to one of your computer's USB ports.

1




Extract our program (Right click → Extract all) then launch « nano\_phyphox\_v1.ino ».

2



3

Go to Tools → Board → Boards Manager, then use the search bar to find and install the last version of Arduino Mbed OS Nano Boards by Arduino.



4

Go to Tools → Manage Libraries, then use the search bar to find and install the last versions of these 5 libraries:

- phyphox BLE by RWTH Aachen University
- Arduino\_LSM9DS1 by Arduino
- Arduino\_LPS22HB by Arduino
- Arduino\_HTS221 by Arduino
- Arduino\_APDS9960 by Arduino


Go to Tools → Port and choose the port COMX to which an Arduino Nano 33 BLE is connected.

```
#include <phyphoxBle.h>
#include <Arduino_LSM9DS1.h>
#include <Arduino_LPS22HB.h>
#include <Arduino_HTS221.h>
#include <Arduino_APDS9960.h>

char board_name[] = "Phyphox1"; // to change the name displayed
float choice = 0.0;
float accx, accy, accz, acc;
```

If needed, change the name displayed by the board (mostly useful if there are several boards in the same room).

5



6

Click on Upload ("→" icon under File and Edit).


Wait for the end of the loading, which will be indicated by "Done uploading" at the bottom of the window.



Download the phyphox app by RWTH Aachen on Google Play or App-Store.

Launch the phyphox app, click on "+" then "Add an experiment with a QR code".  
Flash the following QR code.  
(Alternatively, you can go to our site and click on "Phyphox experiments" under Tutorials and program.)

7



8

Select Save all.

In the list of phyphox experiments, ours are now available under "Arduino Nano 33 BLE Sense".


<b>Acc</b>	Acceleration Plots acceleration measured by the Arduino...
<b>Pre</b>	Analog Input Plots A0 analog input.
<b>Gyr</b>	Angular Velocity (Gyroscope) Plots angular velocity measured by the Arduino...
<b>III</b>	Illuminance Plots illuminance measured by the Arduino...
<b>Mag</b>	Magnetic Field Plots magnetic field measured by the Arduino...
<b>Pre</b>	Pressure Plots pressure measured by the Arduino c...
<b>Tem</b>	Temperature and Humidity Plots temperature and humidity measured...

To start using the experiments, power your board, and activate Bluetooth and GPS on your smartphone.


Launch phyphox and select an experiment under Arduino Nano 33 BLE Sense.

9

10



Select the board to which the smartphone should connect (Phyphox1 here).



Run the experiment.